


# Artificial Iris

## Drawing



## Technical Information



Model	Artificial Iris	Artificial Iris - Fibre Free
Application	Treatment of Full or Partial Aniridia, Fixation by Scleral Sutures	Treatment of Full Aniridia for Sutureless Implantation
Material	Pigmented Silicone Elastomer <b>with</b> Polymer Fibre Meshwork	Pigmented Silicone Elastomer <b>without</b> Polymer Fibre Meshwork
Design	Design Based on the Photo Print-Out (Hard Copy) of the Patient's Eye Serving as Colour Target for the Individual Production with Structure and Topography Resembling the Natural Iris.	
Availability	Delivery Time 4 to 8 Weeks after Receipt of Prescription and Approved Photography.	

## Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13  
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0  
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de  
www.dr-schmidt-iol.de